



1950nm SLD Benchtop Light Source



- **Product Description**

Idealphotonics' modular laser control is based on an advanced microprocessor control system, combined with high-precision ATC and ACC (APC) control circuits to achieve highly stable laser output, while ensuring that the light source is quick and intuitive to operate. We can also provide corresponding communication interfaces and control software according to user requirements to achieve computer control. This light source uses a one-key recovery function (Run/Stop button) to effectively help customers return to the previous working state. This is a highly integrated modular

system light source that uses PC-side software intelligent control. Customers can set the required working temperature and current according to their needs. It is very suitable for experimental scientific research and production testing. In addition, we need to modulate the laser in some application fields. We have connected two external modulation ports, one for high frequency and one for low frequency, to better meet customers' needs for multiple uses of one machine.

- **Product features**

Support one-key restore function (no need to restart and preheat) 、 Software remote control, intelligent control 、 Stable output power, continuously adjustable、 Compact structure、 High-precision ACC and ATC control circuit、 Built-in high and low modulation bandwidth BNC interface

- **Part Number**

MP-SLDS-1950-5-60-SA-M

- **Application area**

Fiber optic transmission system、 Fiber optic gyroscope、 Fiber optic sensor、 Optical coherence tomography、 Test light source



● Core parameters

Wavelength	Spectral Width	Output Power
1950nm	60nm	5mW

● General Parameters

Driver Parameters

characteristic	Min.	Max.	Unit	Note
Power supply voltage	4.8	6.1	VAC	DC, it is recommended to use a regulated power supply
Power	0	10	W	
Laser driver current	0	149 378 624	mA	Patch cord optional
Current resolution (LSB)	Max current/65536		V	
Response frequency	0	25	MHz	-3db



Temperature control range	-10	50	°C	
Temperature resolution	0.001		°C	
Temperature stability		0.002	°C	Requires enclosure cooling
TEC input current	-1.5	1.5	A	
TEC input voltage	-4.4	+4.4	V	
Analog input	0	2.5	V	

Laser parameters

Parameters	Symbol	Min.	Typical value	Max.	Unit
Wavelength	λ		1950		nm
Spectral width	$\Delta\lambda$	50	60	70	nm
Threshold current	I_{th}		30	40	mA
Operating current	I_{op}		200	300	mA
Output power	Pf	3	5	10	mW



Fiber type	SM1950				
Operating voltage	V_f		1.8	2.5	V
Thermistor	RT	9.5	10	10.5	KΩ
Connector type	FC/APC				


Laboratory measured photoelectric data

Optical and Electrical Data						
No.	Parameter	Symbol	Unit	Testing result	Requirement	Condition
1	Threshold current	I_{th}	mA	64	≤80	at 25.0°C
2	Bias during operation	I_o	mA	230	≤230	at 4.65mW power
3	Bias voltage	V_f	V	1.178	≤3	at 4.65mW power
4	Slope efficiency	DE	mW/mA	0.03	0.01-0.05	at 4.65mW power
5	Operating temperature	T	°C	25	15-35	at 4.65mW power




6	Thermistor resist	R	Ohm	10000	9500~10 500	at 25.0°C
7	wavelength	λ	nm	1954.80	1950 ± 10	at 4.65mW power
NOTE: I_{max}=230mA						
Laser Diode characteristic graph				Optical mapping and data		

Test spectrum



Power curve




The burr is water peak absorption, and the actual emission spectrum is smooth






Control software interface




Beam quality





Pin Definition



1	Thermoelectric Cooler(+)
2	Thermistor
3	NC
4	NC
5	Thermistor
6	NC
7	NC
8	NC
9	Laser Cathode(-)
10	Laser Anode (+)
11	NC
12	NC
13	Case Ground
14	Thermoelectric Cooler(-)