




## 2um Single Frequency Single Mode Thulium Doped Fiber Amplifier (TDFA) 1900-2050nm



### ● Product Description

Idealphotonics Thulium-doped Benchtop Single-frequency Single-mode Fiber Amplifier is a power amplifier specially designed for ultra-narrow linewidth, single-frequency signal light (such as fiber lasers based on DFB and DBR principles). This amplifier can amplify low-power optical signals in the kHz range to an output power of up to 50W, and can well maintain the

spectral characteristics of the input signal light. Our amplifier uses a high-power, high-performance multi-mode pump source and double-clad fiber amplification technology, and the output power is continuously adjustable. The benchtop fiber amplifier is a complete Turn-Key system with internal microprocessor control. The front panel is equipped with a laser start switch and a power status LCD display, as well as an output power adjustment knob. Benchtop single-mode fiber amplifier (TDFA) can be widely used in scientific research, coherent synthesis, coherent detection sensor systems, etc.

## ● Product features

High output power、 Low noise figure、 Turn-Key system、 ACC/ALC mode operation

## ● Part Number

MP-TDFA-2000-1W-SM

## ● Application area

Optical communication、 Test and measurement、 Scientific applications

## ● Core parameters

Wavelength	Output Power	Connector
1900-2050nm	1W	FC/APC

## ● General Parameters

### Parameter

PN#	MP-TDFA-2000-1W-SM
	2um fiber amplifier (TDFA)
Application	Pre-amplified optical amplifier
Signal wavelength range	1900-2050nm
Maximum output power*1	>+30dBm />1W
Gain*2	>20dB
Noise factor*2	<6dB
Amplification control	ACC / ALC
Monitoring items	Input power, output power, LD current, LD temperature, ambient temperature
Safety functions	Remote shutdown interlock
Remote control interface	RS232C/IEEE488.2(GP-IB)



Input/output optical fiber	PMF(Nufern PM1950)
Optical fiber connector	FC/APC
Minimum extinction ratio	20dB
Input signal line width	Down to 0.1KHZ
Operating temperature	0 ~ 40 °C
Storage temperature	-10 ~ 60 °C
Dimensions*3 (mm)	88 x 430 x 450
Weight	5 Kg
Power consumption	<30W
Power supply	AC 100 ~ 240V (50/60Hz)

**Notes:**

\* 1 Input power 1-10dBm@1950


\* 2 Under the premise of ensuring beam quality  $M2 < 1.05$ , single mode, TEM00 mode

\* 3 Excluding protruding parts


\* MP-TDFA-2000-1W-SM is designed for 1950~2050nm, total power 0dBm input signal to achieve the best gain stability.



### Spectrum



### Power stability curve and modulation linearity test curve





<b>The following was filled out by AOL Labs staff:</b>			
<b>Test date</b>	<b>2021.11.9</b>	<b>Test engineer</b>	<b>Dai Jiahao</b>
<b>Test product</b>	<b>2um continuous fiber amplifier T DFA (XZ2010T264)</b>	<b>Reviewer</b>	<b>Wang Xiuxiang</b>
<b>Purpose of the experiment</b>	<b>Power stability test, spectrum diagram</b>		
<b>Experimental instruments</b>	<b>Power meter, spectrometer</b>		

—: 2000 nm


### 1. Enlarge the noise test chart




In the above figure, Trace A (yellow) is the seed light spectrum, and Trace B (purple) is the spectrum after amplification (after passing through a 2um attenuator), and the same below. The amplified noise measured at 2000nm is 11.343dB.



## 2. Curve of power



## 3. Power stability



1950 nm


## 1. Amplify the noise test chart




The amplified noise measured at 1950nm is 13.705dB.



## 2. Curve of power



## 3. Power stability



λ: 1910 nm


### 1. Amplify the noise test chart



The amplified noise measured at 1910nm is 13.262dB.



## 2. Curve of power



## 3. Power stability

